

# Object independent MRL analysis GUI

- ① Read Dicom into niftiView datasets  
n.dcm  $\rightarrow$  1.nii.gz  
with/without contours.
- ② Re-orient MRL images
- ③ If required, mask image  
contours, threshold, geometry - mask options.
- ④ Rigid registration to reference image  
(CT, MR)  
 $\rightarrow$  Assess registration
- ⑤ Any cropping to ensure no gray patches  
resulting from rigid result with different  
input image sizes.
- ⑥ Run deformable registration (nifty reg based).  
 $\rightarrow$  have option to change the final  
control point grid spacing depending on  
application
- ⑦ Mask deformation field as may be required  
for analysis : - patient contour  
- subsection of phantom volume  
- ...
- ⑧ Import deformation field data into matlab to:
  - Convert x, y, z voxels into dist from iso
  - Calculate total ~~dist~~ distortion in 3D.
  - Calculate max, mean, min & std dev.